

TE7 Ultrasound System

Transducer Family



C5-1s (Convex)

Applications: Adult Abdominal, Nerve, OB/GYN, Pediatrics, Urology, Vascular
FOV: 101°
Convex radius: 60 mm
Biopsy kit: NGB-022, multi-angle, reusable



L12-3RCs (Linear)

Applications: MSK, Nerve, Pediatrics, Small Parts, Vascular
FOV: 38 mm
Biopsy kit: NGB-043, multi-angle, reusable



C5-2s (Convex)

Applications: Adult Abdominal, Nerve, OB/GYN, Pediatrics, Urology, Vascular
FOV: 76°
Convex radius: 51 mm
Biopsy kit: NGB-015, multi-angle, reusable



L12-4s (Linear)

Applications: MSK, Nerve, Pediatrics, Small Parts, Vascular
FOV: 38 mm
Biopsy kit: NGB-007, multi-angle, reusable



C4-1s (Convex)

Applications: Abdomen, Cardiology, OB/GYN
FOV: 96°
Convex radius: 30 mm
Biopsy kit: NGB-036, multi-angle, reusable; CIVCO 698-019, disposable



L11-3VNs (Linear)

Applications: MSK, Nerve, Small Parts, Vascular
FOV: 38 mm
Biopsy kit: NA
eSpatialNavi™: enabled



C11-3s (Micro-Convex)

Applications: Neonatal Abdominal, Neonatal Cardiac, Neonatal Cranial, Nerve, Neurology, Pediatric Abdominal, Pediatric Cardiac, Vascular
FOV: 100°
Convex radius: 15 mm
Biopsy kit: NGB-018, multi-angle, reusable



L9-3s (Linear)

Applications: MSK, Nerve, Small Parts, Vascular
FOV: 43 mm
Biopsy kit: NGB-034, multi-angle, reusable



L14-6Ns (Linear)

Applications: MSK, Nerve, Pediatrics, Small Parts, Superficial, Vascular
FOV: 38 mm
Biopsy kit: NGB-007, multi-angle, reusable



7L4s (Linear)

Applications: Adult Abdominal, MSK, Neonatal Cranial, Pediatrics, Small Parts, Vascular
FOV: 38 mm
Biopsy kit: NGB-007, multi-angle, reusable

**L7-3s (Linear)**

Applications: Adult Abdominal, MSK, Nerve, Pediatrics, Small Parts, Vascular
FOV: 38 mm
Biopsy kit: NGB-007, multi-angle, reusable

**L20-5s (Linear)**

Applications: MSK, Nerve, Small Parts, Superficial, Vascular
FOV: 28.5 mm
Biopsy kit: NA

**L14-6s (Linear)**

Applications: Pediatrics, MSK, Neonatal Cranial, Nerve, Small Parts, Superficial, Vascular
FOV: 25 mm
Biopsy kit: NGB-016, multi-angle, reusable

**L14-5sp (Linear)**

Applications: Intraoperative, MSK, Nerve, Superficial, Vascular
FOV: 25.3 mm
Biopsy kit: NA

**P10-4s (Sector)**

Applications: Neonatal Abdominal, Neonatal Cardiac, Neonatal Cranial
FOV: 90°
Biopsy kit: NA

**P7-3s (Sector)**

Applications: Abdominal, Cardiac, Cephalic, Pediatric
FOV: 90°
Biopsy kit: NA

**SP5-1s (Sector)**

Applications: Adult Abdominal, Adult Cardiac, Adult Cranial, Pediatric Cardiac
FOV: 90°
Biopsy kit: NGB-011, multi-angle, reusable

**P4-2s (Sector)**

Applications: Adult Abdominal, Adult Cardiac, Adult Cranial, Pediatric Cardiac, Pediatrics, Neonatal Cranial
FOV: 90°
Biopsy kit: NGB-011, multi-angle, reusable

**V11-3ws (Endocavity)**

Applications: OB/GYN, Urology
FOV: 173°
Convex radius: 11 mm
Biopsy kit: NGB-004, single-angle, reusable

**P7-3Ts (TEE)**

Applications: Adult Transesophageal Echo
FOV: 90°
Biopsy kit: NA

**7LT4s (Linear)**

Applications: Adult Abdominal, Intraoperative, MSK, Neonatal Cranial, Pediatrics, Small Parts, Superficial, Vascular
FOV: 43 mm
Biopsy kit: NGB-010, multi-angle, reusable

Mindray North America Innovation Center

2100 Gold Street
 San Jose, CA 95002

Tel: 800.288.2121 Support: 877.913.9663 www.mindray.com

©2018 Mindray DS USA, Inc. Subject to change. Mindray® is a registered trademark of Shenzhen Mindray Bio-Medical Electronics Co. Ltd. 10/18 P/N: 0002-08-40352 Rev B

mindray
 healthcare within reach