

# ME8 Ultrasound System

## Transducer Family



### SC5-1Ns (Convex)

**Applications:** Abdomen, Gynecology, Obstetrics, Vascular, Nerve, Musculoskeletal, Urology, Thoracic/Pleural, Small Parts  
**FOV:** 73°  
**Convex radius:** 60 mm  
**Needle guide:** NGB-022, stainless steel, detachable



### SP5-1Ns (Phased)

**Applications:** Abdomen, Gynecology, Obstetrics, Cardiac, Pediatric, Vascular, Thoracic/Pleural, Cephalic  
**FOV (Max):** 90°  
**Needle Guide:** NGB-011, stainless steel, undetachable



### L12-3RCs (Linear)

**Applications:** MSK, Nerve, Pediatrics, Small Parts, Vascular  
**FOV:** 38 mm  
**Needle guide:** NGB-043, stainless steel, detachable; NGB-044, stainless steel, detachable, out-of-plane



### V11-3Hs (Endocavity)

**Applications:** Gynecology, Obstetrics, Urology  
**FOV:** 170°  
**Needle guide:** NGB-025, stainless steel, undetachable



### L11-3VNs (Linear)

**Applications:** MSK, Nerve, Small Parts, Vascular  
**FOV:** 38 mm  
**Biopsy kit:** NA eSpatialNavi™: enabled



### L20-5s (Linear)

**Applications:** Musculoskeletal, Nerve, Small Parts, Superficial  
**FOV:** 29 mm  
**Biopsy kit:** NA



### L14-6Ns (Linear)

**Applications:** MSK, Nerve, Pediatrics, Small Parts, Superficial, Vascular  
**FOV:** 38 mm  
**Needle Guide:** NGB-007, plastic, detachable

#### Mindray North America

800 MacArthur Blvd., Mahwah, NJ 07430

Tel: 800.288.2121 Support: 877.913.9663 Fax: 800.926.4275 [www.mindray.com](http://www.mindray.com)

Mindray® is a registered trademark of Shenzhen Mindray Bio-Medical Electronics Co., Ltd. All brands and product names are trademarks of their respective owners. ©2020 Mindray DS USA, Inc. Subject to change. 09/20 P/N: 0002-08-40468 Rev A

**mindray**  
NORTH AMERICA